



# PRESSURE INFUSION BAGS

MESH PANEL | CLEAR PANEL | REUSABLE CLEAR PANEL



SINGLE PATIENT USE BAGS

REUSABLE BAGS



**FIGURE 1:** Solution bag hook enables easy bag changing without having to remove the infusion bag from the IV pole; keeps fluid bag in proper position while in use



**FIGURE 2:** Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Clear panel bags enables full visibility of solution bag to monitor fluid level and easily scan patient barcodes



**FIGURE 3:** 360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

8807, 8808, 8809 MESH PANEL  
8810, 8811 CLEAR PANEL  
8817, 8818, 8819 REUSABLE CLEAR PANEL

## QUICK REFERENCE GUIDE

1. Remove pressure infusion bag from package, unfold pressure bag and hang pressure bag on to IV pole by the loop located at the top of the pressure bag (Figure 1).
2. Insert intravenous bag between mesh panel and pressure bag on 8807, 8808 and 8809 only.
3. If using 8810, 8811, 8817, 8818 or 8819, place intravenous bag between clear panel and pressure bag.
4. Hang loop of intravenous bag onto welded solution bag hook located just above panel on pressure bag (Figure 1).
5. FOR 8809 ONLY: Slide hanger loop through top of 3000 ml intravenous bag and rehang on IV pole.
6. Inflate pressure infusion bag by rotating the off portion of the valve/stopcock to the right. Pump inflation bulb until gauge reaches the proper pressure (300 mmHg) indicated in the green portion on the manometer (Figure 2).
7. Once proper pressure is reached on gauge, rotate off portion of the valve/stopcock straight up towards the pressure bag to maintain pressure in bag. **NOTE: Precise positioning of valve/stopcock is important or bag could lose pressure.**
8. To deflate bag turn the off portion on the valve/stopcock straight down towards inflation bulb to release pressure in bag.
9. Proper valve positioning is illustrated on all pressure bags.

**CAUTION: All PremierPro™ Pressure Infusion bags come equipped with a 330 mmHg pressure relief valve that prevents over inflation. Once 330 mmHg is reached, the manometer will indicate red on the manometer along with the word “caution” and pressure will begin to release from valve. Pressure will continue to release from valve until pressure decreases below 330 mmHg (Figure 3).**