

# PremierPro<sup>™</sup> MEMBER SURGICAL GLOVE BROCHURE

Contract Number: PP-OR-2039/AD-OR-2039

www.ansell.com





### 05 About Ansell

- Emergency Preparedness
- Sustainability

### 08 Why Ansell

- Ansell Cytostatic Permeation Program (ACPP)
- Acceptable Quality Level (AQL)
- Worker Experience Innovation (WEI Technologies)
- AnsellCARES<sup>SM</sup>

### **14 Surgical Gloves Overview**

- Material Guide

### **17 Surgical Gloves Portfolio**

- Polyisoprene
- Specialty
- Neoprene
- Latex

### **28 Resources**

- Glove Size Finder
- Why and How to Double Glove

### **ABOUT ANSELL**

Ansell has been protecting people for over 130 years. We have created specialized teams to focus on the personal protection needs of workers in many industries.

# We help patients and healthcare workers stay safe

As the healthcare industry evolves, the pressure to improve patient outcomes and drive efficiencies continues to rise. Ansell, a market-leading innovator in high-quality protection solutions for workplace safety, is committed to help you stay safe and improve outcomes for both patients and healthcare workers.

We are passionate about protecting people, setting ourselves ambitious goals for transformation and to protect sustainably both employees and the environment we live in.

Our global presence ensures we're able to meet the ever-growing demand for surgical gloves. We have 4 locations throughout Asia that are geared to producing the surgical gloves that you trust. Our facilities are paired with a sterilization facility to ensure that every pair of gloves that is donned by surgical professionals are clear from contaminants.

### **Ansell Emergency Preparedness**

A top priority for us is emergency preparedness. We ensure that production, energy, raw materials, and quality assurance are always available.

We place great emphasis on assurance, factory and safety stock, and customer service level to ensure that we can respond quickly to customer needs, providing them with the solutions they require, when they need them.



### **OUR SURGICAL FOOTPRINT**



γ	<b>3</b> IN-HOUSE GAMMA STERILIZATION
$\left( \right)$	Kovai
(	Sri Lanka
(	Malaysia



Diversification in Manufacturing Facilities and Product Lines




Raw Material Multi-Sourcing Strategy



Robust Sales & Operation Planning

### SUSTAINABILITY GOALS & INITIATIVES

At Ansell, we've set ourselves ambitious goals to transform our company in order that we can better protect our employees and the environment we live in.



### Kovai

- \$80 Million investment in the first and foremost Indian greenfield project
- The facility is designed according to Indian Green Building Council (IGBC) criteria
- Built with low-carbon, low-impact, non-toxic materials, and re-purposed materials



### Biyagama

- 83% of the facility runs on renewable energy from biomass<sup>1</sup>
- Facility certified by Intertek for diverting >99% of waste from landfill<sup>2</sup>

Melaka

landfill<sup>2</sup>

• Facility runs on 92% green


Electricity Tariff Program<sup>1</sup>

electricity from onsite rooftop

solar panels and a local Green

• Facility certified by Intertek for

diverting >99% of waste from



<sup>1</sup> 2023 Ansell Sustainability Report; June 2023, ISO 14064-1: 2018 external audit by Control Union, ref.period 7/2022 - 6/2023 ertek, Zero Waste to Landfill certificate

### Packaging Waste: SMART Pack<sup>™</sup> 2.0

#### Ansell has improved its packaging material over the last couple of years to result in a **34%** reduction of CO<sub>2</sub>e per pack





- Smaller pouches, smaller glove boxes, smaller cartons, This all scales down to less plastic and cardboard waste, as well as a **reduced carbon footprint** thanks to lower CO<sub>2</sub> emissions.
- Ansell's more compact dispenser boxes are between 23% to 80% smaller in square area compared to standard packaging in the market.
- Dispenser boxes and cartons are made from 70% and 90%. recycled materials respectively All parts of surgical packaging



\* Always check your local recyclable status of HDPE and LDPE as these plastic materials may not be considered suitable for recycling in your country.

### **Collaborations**

now recyclable.

### °CLIMATE GROUP

Ansell aligns with **Climate Group's** leading global initiatives for a low-carbon future, RE100 and EP100, making good on our climate commitments announced early this fiscal year.



Ansell joined the Renewable Thermal Collaborative, one of the largest coalitions addressing the demand for decarbonization of thermal energy generation.



The Responsible Glove Alliance (RGA) was launched, with seven founding members, including Ansell. The RGA is supported by the **Responsible Business Alliance (RBA)**, the world's largest industry coalition dedicated to corporate social responsibility in global supply chains.

### People

Ansell condemns all forms of modern slavery and human rights abuses, including the use of forced labor. We consistently communicate to suppliers our expectations for having effective systems in place for ensuring the highest standards of health and safety of all workers and remediating human rights abuses identified in their supply chains. To date, we have seen improved labor standards amongst our suppliers.

Ansell recognizes that modern slavery is a profoundly complex issue that requires collaboration and time to rectify and effect lasting industry-wide changes. Therefore, we remain committed to working with the entire international community, industry partners, governments, and non-government organizations to promote a coordinated approach to addressing the root causes of human rights abuses in our industry.

## 

Front:

HDPE plastic

(white)



Back:

LDPE plastic

Surgical poly-pouches will be made using a combination of high density polyethylene (HDPE) and low density polyethylene (LDPE) thermoplastic. HDPE and LDPE are know for their high recyclability, durability, as well as their weather and chemical resistance.

### **Taking Our Information Online**

- Reducing printed materials by giving customers a one stop access point to all product information online.
- Removal of collecting any manual feedback from customers. Customers can now scan a OR code to provide feedback or contact Ansell 24/7.
- Ensuring customers have direct access to Ansell 24/7.



### **Awards & Certifications**



## SILVER 202 ecovadis

### EcoVadis Silver Medal, Stamp of **Excellence in Sustainability**

Ansell was awarded a silver medal by EcoVadis based on reporting data for FY21, placing it in **top 25%** of all companies assessed and the **top 6%** of manufacturers within our applicable industry.



#### WSH Asia Award 2022: **Best Product/Service**

SMART Pack<sup>™</sup> received an excellence recognition within the safety industry based on new and innovative solutions that improve workplace health and safety.

## ANSELL CYTOSTATIC PERMEATION PROGRAM

### **Dynamic Permeation Testing for Greater Confidence at Work**


Ansell's Cytostatic Permeation Program (ACPP), a unique dynamic permeation test, has been designed to give you an added in-use perspective of permeation detection to help you select the right glove for safer handling of chemotherapy drugs.

While ASTM and EN standards determine glove permeation under static test conditions, ACPP mimics the everyday use in your workplace through controlled dynamic conditions.

Current standards testing does not consider all working conditions that may influence the permeation of your glove's protective barrier. Such influencing factors are for example the concentration and exposure time of the chemotherapy drug being handled, glove properties such as thickness; the stretching and flexing motion of the task, as well as both body and workplace temperatures.<sup>4</sup>

It is important you consider how static (standards) vs dynamic (ACPP) testing differs to gain a comprehensive view of how the permeability of a glove to a chemotherapy drug is affected under the different conditions.

### ACPP mimics in-use conditions, and together with Standards, ensures the right glove for safer handling of hazardous drugs.



#### References

- 1. EN stated in this graphic refers to EN 16523-1.
- ASTM stated in this graphic refers to ASTM D6978
- ACPP uses the ADTD unique to this program. For 400+ mm length gloves, both palm and cuff must be tested.
- 5. Detection limit is based on the permeation rate expressed as the amount (in µg) of chemical per surface area (cm²) of the test specimen per minute (min).
- 6. ACPP uses highly sensitive analytical methods, Liquid Chromatography Mass Spectrometry (LS-MS) and Inductively Coupled Plasma Mass Spectrometry (ICP-MS), to enable permeation to be detected at a very low level.
- 7. HIPEC = Hyperthermic intraperitoneal chemotherapy; a highly concentrated, heated chemotherapy treatment
- 8. It is not mandatory to test using a chemotherapy drug for EN 16523-1.
- All 9 drugs tested for ASTM D6978 are chemotherapy drugs.
- 10. All 27 drugs tested for ACPP are chemotherapy drugs.

## **ACCEPTABLE QUALITY LEVEL**

Quality medical gloves are essential to protect healthcare workers and patients from potential hazards, and can be measured by Acceptable Quality Level.

### What is Acceptable Quality Level (AQL)?

AQL is an industry standard that is a statistical sampling process for evaluating quality. According to the International Standards Organization (ISO) (2859-1: 1999)<sup>1</sup>, AQL is "the worst tolerable process average when a continuing series of lots is submitted for acceptance sampling". Process average is the typical percentage of defective gloves in the lots/batches sampled.

### What are the required AQLs?

Various international standards determine the AQL that manufacturers must comply to. However, manufacturers can set their own standards as long as they are stricter than the international standards. Ansell has these AQL Standards:

- Ansell surgical gloves adhere to an AQL of 0.65 (Meeting or exceeding recommendations set by ASTM, EN & ISO)
- The lower the AQL, the higher the level of glove barrier protection and the lower the incidence of finding a defect resulting in a higher quality product.
- Technical data sheets containing the AQL for each specific glove are available upon request.

### How is AQL determined and measured for gloves?



#### Water Leak Test to determine the AOL level of a glove for pinholes

- Gloves are filled with 1000 ml (1 litre) of water.
- They are bound or sealed and hung upside down for two minutes
- Leaks are checked while gloves are under sustained pressure.
- This is the recognized test method for global glove standards.

### Ansell's Acceptable Quality Level, a guarantee for higher quality and protection.

Medical gloves are critical to protect both patients and healthcare workers and as such, they must be manufactured to a very high standard at all times. Amongst the most critical protective features are the absence of pinholes. The accepted rate of defects allowable to be released after passing though control processes is measured by AQL, and therefore this test is one of the most direct indicators of manufacturers' quality and process standards.

#### References:

- http://www.iso.ora/iso/home/store/catalogue\_tc/catalogue\_detail.htm?csnumber=1141 Accessed 26/11/2015.
- 3. American Society for Testing and Materials, ASTM D3577-19 Standard Specification for Rubber Surgical Gloves http://www.astm.org/Standards/D3577.htm Accessed 08-03-22.

- 6. Japanese Industrial Standard, JIS T 9107:2018 Single-use Sterile Surgical Rubber Gloves <a href="http://infostore.saiglobal.com/store/Details.aspx?ProductID=800132">http://infostore.saiglobal.com/store/Details.aspx?ProductID=800132</a> Accessed 4-12-15.

#### Table 1: **International Surgical Glove Standards**

Surgical Glove Standards	Inspection Level	AQL
AS/NZS 4179: 2014 <sup>2</sup> Applicable to Australia/ New Zealand	G1	1.0
ASTM D3577-19 <sup>3</sup> Applicable to US & Canada	G1	1.5
EN 455-1:2020 <sup>4</sup> Applicable to the European Union	G1	0.65
ISO 10282: 2014 <sup>5</sup> Adopted by the rest of the World	G1	1.5
JIS T9107: 2018 <sup>6</sup> Applicable to Japan	G1	1.5

AQL is a pass/fail where a predetermined sample size of a manufactured lot is tested following the sampling plan and protocols established by the various international standards or more stringent standards set by manufacturers to ensure stricter and higher quality is delivered to the customer.

The sampling plan is an inspection procedure of a sample size (which are gloves that are randomly selected from a batch of gloves) that is used to determine acceptance or rejection criteria from an inspection batch or lot. The sample size to be tested is set by:

#### The lot or batch size

The inspection level - determined by region as reflected in Table 1

The decision on whether or not to release the full production batch to the market will be based on the results obtained from the sample size tested, as set by the AQL value specified in standards for the market, or determined by individual manufacturers. Manufacturers will also decide whether to discard, or rework the lot until the AOL is achieved.

Standard New Zealand Te Mana Tautikanga O Aotearoa, AS/NZS 4179:2014 Single-use sterile rubber surgical gloves http://infostore.saiglobal.com/store/Details.aspx?ProductID=1742680 Accessed 4-12-15 4. European norm, EN 455 Part 1: 2020 Medical Gloves for single use http://ec.europa.eu/growth/single-market/european-standards/harmonised-standards/ medical-devices/index\_en.htm Accessed 08-03-22. 5. International Standardization Organization, ISO 10282:2014 Single-use sterile rubber surgical gloves http://www.iso.org/iso/home/store/catalogue\_tc/catalogue\_te/catalogu

<sup>1.</sup> International Standards Organization, ISO 2859-1:1999 Sampling procedures for inspection by attributes — Part 1: Sampling schemes indexed by acceptance quality limit (AQL) for lot-by-lot inspection

### **WEI TECHNOLOGIES**

### Through our Worker Experience Innovation (WEI) program, we develop a deeper understanding of workers' day-to-day experiences and challenges. We then transform these insights into new WEI technologies and product ideas.



### **ANSELL PROPRIETARY TECHNOLOGIES**

#### **{>** A.R.T Pinhole Barrier Technology

This technology is proprietary process, integrated into the manufacturing line that helps add an extra level of robustness by minimizing possible defects that may later result in pinholes.



A coating is applied over the glove former to seal any micro-defects, creating a smooth finish and delivering a consistent thickness across the glove. This minimizes irregularities that may later result in pinholes

#### **GLOVE STYLES**

- Ø PremierPro<sup>™</sup> Polyisoprene
- Ø PremierPro<sup>™</sup> Polyisoprene Underglove
- Ø PremierPro<sup>™</sup> Polyisoprene Micro
- Ø PremierPro<sup>™</sup> Polyisoprene Hybrid
- Ø PremierPro<sup>™</sup> Polyisoprene Ortho
- Ø PremierPro<sup>™</sup> Neoprene Micro
- Ø PremierPro<sup>™</sup> Latex
- 🥑 PremierPro<sup>™</sup> Latex Ortho
- ⑦ PremierPro<sup>™</sup> Latex Underglove
- ⑦ PremierPro<sup>™</sup> Latex Micro
- GAMMEX<sup>®</sup> Non-Latex PI Textured
- 🧿 GAMMEX® PI Plus Glove-in-Glove<sup>™</sup> System


PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039

- GAMMEX<sup>®</sup> Non-Latex
- 6 ENCORE<sup>®</sup> Latex Textured

#### DERMASHIELD<sup>®</sup> Donning Aid Technology

This technology is a proprietary inner coating for fast and easy damp or dry hand donning whether it is the first don or during intraoperative changes.





Hydrophilic lubricants in the coating are activated upon contact with wet or damp hands to facilitate damp donning and help minimize the feel of sweaty hands

#### **GLOVE STYLES**

- Ø PremierPro<sup>™</sup> Polyisoprene
- ⑦ PremierPro<sup>™</sup> Polyisoprene Underglove
- Ø PremierPro<sup>™</sup> Polyisoprene Micro
- Ø PremierPro<sup>™</sup> Polyisoprene Hybrid
- ⑦ PremierPro<sup>™</sup> Polyisoprene Ortho
- Ø PremierPro<sup>™</sup> Neoprene Micro
- Ø PremierPro<sup>™</sup> Latex
- Ø PremierPro<sup>™</sup> Latex Micro
- ⑦ PremierPro<sup>™</sup> Latex Ortho
- Ø PremierPro<sup>™</sup> Latex Underglove
- **GAMMEX**<sup>®</sup> Non-Latex PI Textured
- 🧿 GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System
- GAMMEX<sup>®</sup> Non-Latex

#### ENCORE<sup>®</sup> Latex Textured

#### This technology enables Ansell to design safety solutions that support musculoskeletal health during repetitive tasks to improve worker performance.

Ergonomic Design Technology

ERGOFORM


Supports musculoskeletal health by reducing muscle effort and strain on joints, ligaments, and tendons




Ergonomic design enhances flexibility for maximum range of motion and superior comfort

#### **GLOVE STYLES**

- Ø PremierPro<sup>™</sup> Polyisoprene
- Ø PremierPro<sup>™</sup> Polyisoprene Underglove
- Ø PremierPro<sup>™</sup> Polyisoprene Micro
- Ø PremierPro<sup>™</sup> Polyisoprene Hybrid
- Ø PremierPro<sup>™</sup> Neoprene Micro
- Ø PremierPro<sup>™</sup> Latex
- Ø GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System
- ENCORE<sup>®</sup> Latex Textured





Traditional PI Products

Standard chemical accelerators may not remain within the glove film and its residue, upon contact with skin, could result in Type IV chemical allergies and sensitivities.

#### **GLOVE STYLES**

⑦ PremierPro<sup>™</sup> Polyisoprene

- ⑦ PremierPro<sup>™</sup> Polyisoprene Ortho
- GAMMEX<sup>®</sup> Non-Latex PI Textured

### MANUFACTURING PROCESS TECHNOLOGIES

waste and space -DERMASHIELD<sup>™</sup> Technology

### **GLOVE STYLE**

### **ANSELL PROPRIETARY TECHNOLOGIES**

⑦ PremierPro<sup>™</sup> Polyisoprene Underglove

- Ø PremierPro<sup>™</sup> Polyisoprene Micro
- Ø GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System

are used.

#### **GLOVE-IN-GLOVE**<sup>™</sup> Pre-donned Double Gloving System

#### GAMMEX<sup>®</sup> PI Plus Glove-In-Glove™ System

World's first-to-market innovative pre-donned double gloving system. Simplified process increases speed of donning, and encourages double gloving



Ø GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System

#### SENSOPRENE Ultra Thin Neoprene Technology

This technology delivers a glove up to 30% thinner than Ansell's standard neoprene surgical gloves so wearers' experience a high degree of sensitivity without compromising on durability and the optimum skin protection needed from allergies and sensitivities caused by latex or chemical accelerators.



Ansell's proprietary technology results in a thinner glove without compromising on strength

#### **GLOVE STYLE**

⑦ PremierPro<sup>™</sup> Neoprene Micro

#### SUREFIT Cuff Grip Technology

SUREFIT<sup>™</sup> Technology ensures the inner cuff edge retains its natural surface friction to enable it to securely grip to the surgical gown, helping to eliminate cuff roll-down and heighten protection.



- PremierPro<sup>™</sup> Latex Micro
- PremierPro<sup>™</sup> Latex Underglove
- ⑦ PremierPro<sup>™</sup> Polyisoprene
- Ø PremierPro<sup>™</sup> Latex Underglove
- GAMMEX<sup>®</sup> Non-Latex PI Green GAMMEX<sup>®</sup> Non-Latex PI Textured
- ⑦ GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System
- GAMMEX<sup>®</sup> Non-Latex
- ENCORE<sup>®</sup> Latex Textured



As part of our long-standing commitment to advancing education, AnsellCARES<sup>SM</sup> was created with the mission of adding value beyond product.

Our global multifaceted education programs are guided by clinical experts from around the world. Training resources are based on the latest evidence and practice guidelines to promote safety and keep you informed of the latest developments for application to practice.

Resources are available in several formats to meet your learning needs.

### WHAT WE DO

## EDUCATION

### Industry-leading learning experiences

## Accredited Continuing Education (CE) Programs

Approved provider of accredited CE programs and/or non-accredited courses for a variety of healthcare professionals.

eLearning with HealthStream

(also compatible with other LMS platforms)
Partners in Protection Webinar Series

- (Live and On-Demand)
- Live and Virtual Presentations
- Self-Study Modules



### Grounded in clinical research

Content based on scientific evidence, best practice standards, and new industry innovation trends focused on immediate application to practice.

- Clinical Perspectives
   Knowledge sharing from practicing clinicians
- **Clinical Bulletins** Literature reviews of multiple topics
- **Clinical Study Summaries** Short reviews focused on application to practice
- Knowledge Flash Topics One-minute highlights for busy healthcare professionals

## ENGAGEMENT

### **Connected with healthcare and safety leaders**

Clinical partnerships that offer safety assessment services and multimedia resources designed to promote safety, support risk reduction, and improve outcomes.

#### AnsellGUARDIAN®

On-site consultative O.R room turnover and patient transfer and repositioning safety assessments with recommendations and customized training based on results to help your teams improve safety, clinical outcomes, efficiency and cost

Product Training Resources

- eLearning with HealthStream (also compatible with other LMS platforms)

- Informational Videos
- Step by Step Posters
- InTouch Communications Newsletters. One Page InService Briefs and Blogs



### HOW WE DO IT

### Surgical Glove Evaluation and Implementation\*

Dedicated, O.R. Experienced Clinical Education Managers, Clinical Implementation Teams, and Account Managers

When it comes to conversion support, you can rely on Ansell's team of highly skilled professionals to make this transition a success. Our experts will be with you every step of the way, to assist, guide and support you through this important change process.

Our PROVEN, simple and effective evaluation and implementation process ensures long-term success with minimal disruption for your staff and surgical workflow.



L	
نی	5

### Differentiated Process with 4-Phases that Support a Seamless Transition for Your Team.

**Proven Process** 4 Phases Of Differentiation

#### Assessment

Customization based on the needs of your surgical and procedural teams

### Efficiency

Minimal disruption to staff and workflow procedures

# CONVERSION CONFIDENCE

### Satisfaction

Highly engaged and satisfied staff upon conversion completion

### Follow Up

Available to support your long-term success and provide ongoing educational support

Scan or <u>click here</u>

for more info

PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039

**96%** of conversion participants rated Ansell **superior or acceptable** compared to a competing product\*

97%
 overall satisfaction rate
 with Ansell's evaluation
 and implementation
 process\*

\*Ansell Healthcare, LLC, data on file.

### **Glove Sizing**

Ensuring your team is wearing their preferred style for optimal protection, performance and comfort

#### Engagement

Working alongside staff to collect feedback throughout the entire process

### Optimization

Improving glove quality, outcomes, smart packaging and SKU reduction

### Communication

Open dialogue and partnership with your staff make a successful transition

# SURGICAL GLOVES OVERVIEW

### LOW ALLERGENICITY

Skin-friendly polyisoprene (PI) gloves free from chemical accelerators known to cause chemical allergies, 100% chemical accelerator-free neoprene synthetics and low-protein latex range minimize your risk of developing sensitivities and triggering allergic reactions that can hinder your performance.

### **ERGONOMIC HEALTH**

Exclusive formers ergonomically engineered and certified\* to custom-fit the natural shape of your hands reduce the risk of repetitive stress and muscle fatigue. Innovative, proprietary formulations are softer and thinner yet remarkably strong, optimizing comfort without compromising barrier effectiveness.

\*select products have been certified by U.S. Ergonomics.

### **USER-FOCUSED TECHNOLOGIES**

Proprietary technologies are designed to enrich your user experience and enhance performance—including skin moisturizing, damp donning, secure cuffs, enhanced grip, double gloving systems and more.

### **ECO-FRIENDLY PACKAGING**

Patented, environmentally-friendly SMART Pack<sup>™</sup> design is a unique and compact design that focuses on how our packs work for you in the O.R while having less impact on our planet.




### **MATERIAL GUIDE**

Materials matter, especially in the operating room. The right glove material can make the difference between an effortless procedure and an uncomfortable one. Materials can affect comfort, allergic reactions, and even the duration a glove can be worn. By understanding the properties of glove materials, surgeons can select gloves that provide optimal protection, superior grip, and the utmost comfort, ensuring that their focus remains solely on their life-saving work.

	POLYISOPRENE	NEOPRENE	LATEX
LEVEL OF BARRIER PROTECTION	Excellent	Excellent	Excellent
ALLERGEN CONTENT	Very Good	Very Good to Excellent	Varies
STRENGTH & DURABILITY	Very Good	Excellent	Excellent
ELASTICITY	Excellent	Excellent	Excellent
FIT, COMFORT	Excellent	Very Good	Excellent
ECONOMY	Fair	Good	Very Good

Please note that the following descriptions and tables are representative of basic material characteristics. Individual glove performance may vary.



PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039



# SURGICAL GLOVES PORTFOLIO

### **POLYISOPRENE RANGE**

### **PremierPro<sup>™</sup> Polyisoprene**

Thickness	Standard		
Cuff Style	Beaded with SUREFIT™ Technology		
External Glove Surface	Smooth with micro-textured finish		
Internal Glove Surface	Polymer co DERMASHI	ated with ELD™ Technology	
Double Gloving Recommendation	Outer glove	e or underglove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared		
Freedom From Holes (AQL)	0.65		
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
	SIZE	CODE	
	5.5	43155	
	6	43160	
	6.5	43165	
Order Number	7	43170	
	7.5	43175	
	8	43180	
	8.5	43185	
	9	43190	

### **PremierPro<sup>™</sup> Polyisoprene Underglove**

Thickness	Micro		
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology		
External Glove Surface	Smooth with micro-textured finish		
Internal Glove Surface	Polymer co DERMASHIE	ated with ELD™ Technology	
Double Gloving Recommendation	Underglove		
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)		
Freedom From Holes (AQL)	0.65		
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
	SIZE	CODE	
	5.5	43255	
	6	43260	
Order Number	6.5	43265	
	7	43270	
	7.5	43275	
	8	43280	
	8.5	43285	
	9	43290	

### and hand protection Smooth with micro-textured finish: Excellent instrument handling and grip \*No evidence of sensitization reported in independent Human Repeat Insult Patch Test (HRIPT)

Polyisoprene, skin-friendly surgical glove with a soft

Skin-friendly formulation: Not known to cause allergic

Soft, yet strong: Improved wearer comfort, dexterity

glove formulation delivering extra comfort

### FEATURED TECHNOLOGIES

**KEY FEATURES & BENEFITS** 

contact dermatitis\*

DERMASHIELD<sup>\*\*</sup>



EN ISO 374-1 Type B EN ISO 374-5 ľ \$€ КМРТ VIRUS



Polyisoprene, skin-friendly surgical underglove

Micro thin: Optimal comfort, fit and feel when

Green color: Quick and easy breach detection

Skin-friendly formulation: Not known to cause

\*No evidence of sensitization reported in independent Human Repeat

DERMASHIELD<sup>®</sup>

SUREFIT"

with a high level of tactile sensitivity

allergic contact dermatitis\*

**KEY FEATURES & BENEFITS** 

double gloving

Insult Patch Test (HRIPT)

EN ISO 374-1

Type B

🔶 PI-KARE"

FEATURED TECHNOLOGIES

STANDARDS AND CERTIFICATIONS

EN ISO 374-5

**&** 

VIRUS



ERGOFORM<sup>\*\*</sup>

#### STANDARDS AND CERTIFICATIONS

### PI (Polyisoprene)

Scan or click here

for more info

■₩₩



### Scan or <u>click here</u>

for more info



### **PremierPro<sup>™</sup> Polyisoprene Micro**

ickness	Micro		
ff Style	Beaded with SUREFIT <sup>™</sup> Technology		
ternal Glove rface	Smooth with micro-textured finish		
ernal Glove rface	Polymer C DERMASH	oated with IELD™ Technology	
uble Gloving commendation	Outer glov	e	
sted for Use with emotherapy ugs	No, not tested in accordance with ASTM D6978		
eedom From les (AQL)	0.65		
ckaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
	SIZE	CODE	
	5.5	43555	
	6	43560	
	6.5	43565	
der Number	7	43570	
	7.5	43575	
	8	43580	
	8.5	43585	
	9	43590	



Patch Test (HRIPT)

КМРТ

finesse

### PremierPro<sup>™</sup> Polyisoprene Hybrid

ickness	Standard		Hybrid synt
ff Style	Beaded wit Technology	comfort of	
ternal Glove	Smooth wit	h	KEY FEATUR
rface	micro-textu	red finish	<ul> <li>Hybrid to</li> </ul>
ernal Glove	Polymer Co		durabilit
rface	DERMASHI	ELD <sup>™</sup> Technology	<ul> <li>Low-alle</li> </ul>
ouble Gloving commendation	Outer glove	9	range of
sted for Use with			<ul> <li>Semi-tra detection</li> </ul>
emotherapy ugs	Yes, in acco D6978 and	*Free from natur	
eedom From bles (AQL)	0.65	2-mercaptobenz cetylpyridinium	
ckaging	50 pairs per carton/case carton/case	FEATURED T	
	SIZE	CODE	
	5.5	43355	
	6	43360	HYBRID <sup>™</sup> Synthetic Blend
	6.5	43365	
der Number	7	43370	
	7.5	43375	STANDARDS
			EN ISO 374-1 Type A EN
	8	43380	) <sup>(10</sup>
	8.5	43385	
	9	43390	КМРТ

PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039



PI (Polyisoprene)

## POLYISOPRENE RANGE

#### PI (Polyisoprene)

#### Polyisoprene, skin-friendly surgical glove with exceptional dexterity and tactile sensitivity for surgical

#### **KEY FEATURES & BENEFITS**

• Micro thin: Made 15% thinner\* for superb tactile sensitivity and dexterity

• Non-sensitizing formulation: No known risks of allergic contact dermatitis\*

• Smooth with micro-textured finish: Secure grip of surgical instruments

\*Compared to Ansell's standard polyisoprene (PI) gloves

\*\*No evidence of sensitization reported in independent Human Repeat Insult

#### FEATURED TECHNOLOGIES



**SUREFIT**<sup>®</sup>

DERMASHIELD<sup>®</sup>

VIRU



Scan or click here for more info



#### Polyisoprene/Neoprene

#### /brid synthetic surgical gloves offering the mfort of polyisoprene and strength of neoprene

#### FEATURES & BENEFITS

Hybrid technology: The comfort of polyisoprene and durability of neoprene

Low-allergenic formulation: Reduced risks to a range of allergens\*

Semi-transparent design: Fast and easy breach detection when double gloving

e from natural rubber latex proteins, as well as Zinc ercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and lpyridinium chloride (CPC)



#### ATURED TECHNOLOGIES



#### ANDARDS AND CERTIFICATIONS

EN ISO 374-5

(&

VIRUS





## POLYISOPRENE RANGE

### **PremierPro<sup>™</sup> Polyisoprene Ortho**

Thickness	Ortho		
Cuff Style	Beaded with SUREFIT™ Technology		
External Glove Surface	Smooth wit micro-textu		
Internal Glove Surface	Polymer coa DERMASHIE	ated with ELD™ Technology	
Double Gloving Recommendation	Underglove		
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)		
Freedom From Holes (AQL)	0.65		
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
	SIZE	CODE	
	6	43660	
	6.5	43665	
Order Number	7	43670	
	7.5	43675	
	8	43680	
	8.5	43685	
	9	43690	

#### Polyisoprene, skin-friendly, ultra-thick surgical glove to withstand the rigors of surgery, especially where power tools and bone fragments are present

#### **KEY FEATURES & BENEFITS**

- 50% thicker design\*: Enhanced protection from bone fragments
- Skin-friendly formulation: Reduced risk of allergic contact dermatitis\*
- Micro-textured palm and fingertips: Firm grip of power tools

\*Compared to Ansell's standard powder-free latex gloves

\*\*No evidence of sensitization reported in independent Human Repeat Insult Patch Test (HRIPT)

#### FEATURED TECHNOLOGIES



SUREFIT<sup>®</sup> 🛟 PI-KARE"

allergy protection and superb comfort

dexterity and hand protection

**KEY FEATURES & BENEFITS** 

instrument handling



PI (Polyisoprene)



### PI (Polyisoprene)

for more info



### **GAMMEX®** Non-Latex PI Textured

Thio

Cuf

Exte

Surf

Inte

Surf

Test Che

Dru

Free

Hol


ckness	Standard	Polyisopre	
f Style	Beaded witl Technology	textured fii and grip	
ernal Glove face	Textured fin	KEY FEATUR	
ernal Glove face	Polymer coa DERMASHIE	ated with ELD™ Technology	Textured     or dry co
ted for Use with emotherapy	Yes, in acco D6978 and	Donning     doffing	
gs edom From .es (AQL)	0.4	Skin-frie     dermatit	
kaging	50 pairs per carton/case carton/case	*No evidence of Patch Test (HRIF	
	SIZE	CODE	FEATURED 1
	5.5	46755	🔶 A.R.T
	6	46760	Pinhole Barrie
	6.5	46765	Skin-friendly
ler Number	7	46770	
	7.5	46775	STANDARDS
	8	46780	EN ISO 374-1 Type B El
	8.5	46785	
	9	46790	КМРТ

## GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System

Thickness	Standard		Synthetic, pre-donne			
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology		Beaued with SOREFIT		two polyisoprene (PI to double glove in or	
External Glove Surface	Textured fir	ngers and palm	KEY FEATURES & BENE	FITS		
Internal Glove Surface	Polymer co DERMASHI	ated with ELD™ Technology	Pre-donned double     donning time	e-gloving syst		
Tested for Use with Chemotherapy Drugs	Yes, in acco	rdance with ASTM US FDA cleared	• 50% less inner pack footprint	kaging: Redu		
Freedom From Holes (AQL)	0.65		<ul> <li>Skin-friendly PI form known allergens*</li> </ul>	mulation: No		
Packaging	2 pairs per pouch; 50 pairs per box (25 pouches); 4 boxes per carton/case; 200 pairs per carton/case (100 pouches)		*No evidence of sensitization r Patch Test (HRIPT) FEATURED TECHNOLO			
	SIZE	CODE				
	5.5	46355	4			
	6	46360	Skin-friendly PI Technology	GLOVE-IN Pre-donned Dou		
Order Number	6.5	46365				
Order Number	7	46370				
	7.5	46375	STANDARDS AND CERT	IFICATIONS		
	8	46380	EN ISO 374-1:2016 EN ISO 374-5 Type A			
	8.5	46385				
	9	46390				

## **GAMMEX®** Non-Latex PI Green

Thickness	Standard		
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology		
External Glove Surface	Smooth with finish	n micro-textured	
Internal Glove Surface	Polymer coa DERMASHIE	ated with LD™ Technology	
Double Gloving Recommendation	Outer glove or underglove		
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared		
Freedom From Holes (AQL)	0.65		
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
	SIZE	CODE	
	5.5	46155	
	6	46160	
	6.5	46165	
Order Number	7	46170	
	7.5	46175	
	8	46180	
	8.5	46185	
	9	46190	



cetylpyridinium chloride (CPC)

\*Free from natural rubber latex proteins and do not contain Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and

Soft yet strong formulation: Exceptional comfort,

#### FEATURED TECHNOLOGIES

SUREFIT"

## Scan or click here





PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039

h Test (HRIPT)

A.R.T

```
PI-KARE<sup>®</sup>
```

374-1 pe B EN ISO 374-5 Ň \$ КМРТ VIRUS

Holes (A
Packagir





## POLYISOPRENE RANGE

#### PI (Polyisoprene)

#### visoprene, skin-friendly surgical glove with tured finish for exceptional instrument handling

#### FEATURES & BENEFITS

Textured external finish: Firm instrument grip in wet or dry conditions

Donning aid technology: Fast and easy donning and

Skin-friendly, non-sensitizing formulation: Low dermatitis potential'

evidence of sensitization reported in independent Human Repeat Insult

#### TURED TECHNOLOGIES

DERMASHIELD<sup>®</sup>

SUREFIT<sup>®</sup>

#### NDARDS AND CERTIFICATIONS







PI (Polyisoprene)

#### pre-donned double gloving system of oprene (PI) gloves in one, making it easy glove in one don

#### ES & BENEFITS

ed double-gloving system: Saves valuable

inner packaging: Reduced environmental

ndly PI formulation: No sensitization to llergens\*

ensitization reported in independent Human Repeat Insult



#### **TECHNOLOGIES**



🛟 DERMASHIELD"

GLOVE-IN-GLOVE"

SUREFIT

ERGOFORM<sup>®</sup>

Scan or <u>click here</u> for more info



PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039

### **SPECIALTY RANGE**

### **GAMMEX® PI Radiation Attenuation**

Thickness	Standard	
Cuff Style	Straight	
External Glove Surface	Textured Fingers	5
Internal Glove Surface	Chlorinated	
Double Gloving Recommendation	Underglove	
Tested for Use with Chemotherapy Drugs	No, not tested in accordance with ASTM D6978	
Freedom From Holes (AQL)	1.0	
Packaging	5 pairs per box; 4 boxes per carton/case; 20 pairs per carton/case	
	SIZE C	ODE
	6.5 4	6865
	7 4	6870
Order Number	7.5 4	6875
	8 4	6880
	8.5 4	6885
	9 4	6890

**GAMMEX® Cut-Resistant Glove Liner** 

Standard

N/A

gloves

1.0

Color-coded

Knitted fabric

In between two surgical

with ASTM D6978

5 pairs per box;

5.5

6 - 6.5

7 - 7.5

8 - 8.5

1 box per carton/case

CODE

46950

46960

46970

46980 46990

No, not tested in accordance

Thickness

Cuff Style

Surface

External Glove

Internal Glove Surface

Double Gloving Recommendation

Tested for Use wit

Chemotherapy

Freedom From

Holes (AQL)

Packaging

Order Number

Drugs

#### Polyisoprene sterile gloves, protecting hands from risks of scattered radiation exposure

#### **KEY FEATURES & BENEFITS**

- High attenuation property: 69.1% protection from scattered radiation'
- Optimal thickness: Ideal balance of dexterity and hand protection
- Lead-free, skin-friendly formulation: Eliminates risks of lead and latex

\*At a beam energy level of 60 kVp (data on file)

### STANDARDS AND CERTIFICATIONS





Scan or click here

for more info

Polyisoprene with Tungsten

#### High Density Fiber: Dyneema® Diamond, Nylon, Spandex, Polyester

Non-latex, Dyneema® Diamond yarn cut-resistant glove liners provide adequate cut resistance and protection during rigorous procedures\* **KEY FEATURES & BENEFITS** 

- ANSI A2 cut rating: Cut protection during rigorous, high-risk procedures\*
- Lightweight Dyneema® Diamond yarn: Uncompromising comfort and dexterity
- Allergy optimized formulation: Non-latex and chemical accelerator-free

\*WARNING: No glove liner provides complete protection against cuts. Use with care when handling sharp instruments. The glove liner offers protection against laceration-type and scalpel injuries, but not needle stick . injuries.

STANDARDS AND CERTIFICATIONS



Scan or <u>click here</u> for more info



### **ENCORE®** Radiation Attenuation

Thickness	Standard	
Cuff Style	Straight	
External Glove Surface	Smooth with micro-textured finish	
Internal Glove Surface	Polymer coated	
Double Gloving Recommendation	Underglove	
Tested for Use with Chemotherapy Drugs	No, not tested in accordance with ASTM D6978	
Freedom From Holes (AQL)	1.0	
Packaging	5 pairs per box	
	SIZE	CODE
	6.5	45865
Order Number	7	45870
	7.5	45875
	8	45880
	8.5	45885
	9	45890



### **SPECIALTY RANGE**

Natural Rubber Latex with **Bismuth Oxide** 

#### Latex, powder-free personal protective shield protects hands from scattered radiation exposure during fluoroscopic procedures

#### **KEY FEATURES & BENEFITS**

• Superb attenuation property: 53.7% reduction in skin-dose exposure'

Ideal thickness: Superb balance of hand protection and tactile sensitivity

Micro-textured finish: Better grip control and instrument handling \*At a beam energy level of 60 kVpv







### **NEOPRENE RANGE**

### **PremierPro<sup>™</sup> Neoprene Micro**

Thickness	Micro	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Smooth with finish	micro-textured
Internal Glove Surface	Polymer coa DERMASHIE	ted with LD™ Technology
Double Gloving Recommendation	Outer glove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	5.5	43455
	6	43460
	6.5	43465
Order Number	7	43470
	7.5	43475
	8	43480
	8.5	43485
	9	43490

**GAMMEX®** Non-Latex

Thickness

Cuff Style

Surface

Surface

Drugs Freedom From

External Glove

Internal Glove

Double Gloving

Recommendation

Tested for Use wit

Chemotherapy

Holes (AQL)

Packaging

Order Number

Standard

finish

0.65

carton/case

5.5

6

6.5

7

7.5

8

8.5

9

Technology

Beaded with SUREFIT<sup>™</sup>

Polymer coated with DERMASHIELD<sup>™</sup> Technology

Smooth with micro-textured

Outer glove or underglove

Yes, in accordance with ASTM

D6978 and US FDA cleared

50 pairs per box; 4 boxes per

46455

46460

46465

46470

46475

46480

46485

46490

carton/case; 200 pairs per

### Non-latex, sterile neoprene surgical gloves for greater sensitivity and dexterity

#### **KEY FEATURES & BENEFITS**

- Non-latex, chemical accelerator-free formulation: Optimal allergy protection
- Ultra-thin technology: 30% thinner\* for sensitive yet durable gloves
- Smooth surface with micro-textured finish: Secure instrument grip and handling

\*When compared with standard Ansell neoprene surgical gloves

#### FEATURED TECHNOLOGIES

- 🔶 A.R.T DERMASHIELD<sup>\*</sup>
- SENSOPRENE SUREFIT

### STANDARDS AND CERTIFICATIONS



## Neoprene (Polychloroprene)

Scan or click here

for more info

?- ( **•** 



Scan or <u>click here</u>

for more info

#### FEATURED TECHNOLOGIES

🔶 A.R.T 



### STANDARDS AND CERTIFICATIONS



## Neoprene, chemical accelerator-free surgical

#### **KEY FEATURES & BENEFITS**

Chemical accelerator-free formulation: Advanced protection from allergies

gloves offer advanced allergy protection

- Neoprene film strength: Proven to withstand the most rigorous procedures
- Smooth with micro-textured finish: Confident instrument handling and grip



### **PremierPro<sup>™</sup> Latex**

Thickness	Standard	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Smooth	
Internal Glove Surface	Polymer coa DERMASHIE	ated with ELD™ Technology
Double Gloving Recommendation	Outer glove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTN D6978 and US FDA cleared	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	5.5	44155
	6	44160
Order Number	6.5	44165
	7	44170
	7.5	44175
	8	44180
	8.5	44185

#### STANDARDS AND CERTIFICATIONS EN ISO 374-1 Type B EN ISO 374-5 ľ \$₽ KPT VIRUS

🛟 A.R.T

### **PremierPro<sup>™</sup> Latex Micro**

9

Thickness	Micro	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Smooth wit finish	h micro-textured
Internal Glove Surface	Polymer co DERMASHII	ated with ELD™ Technology
Double Gloving Recommendation	Outer glove or underglove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	5.5	44555
	6	44560
	6.5	44565
Order Number	7	44570
	7.5	44575
	8	44580
	8.5	44585
	9	44590

44190

#### EN ISO 374-1 EN ISO 374-5 Type B ľ ( 🏵

🛟 A.R.T

SUREFIT<sup>\*\*</sup>



## LATEX RANGE

#### Natural Rubber Latex

#### Latex, powder-free surgical gloves for a comfortable. natural fit and feel

#### **KEY FEATURES & BENEFITS**

Soft yet robust formulation: Optimized comfort and tactile sensitivity

Ergonomic design: Reduced hand muscle effort during surgery

Micro-textured external finish: Safe and secure grip of medical instruments



#### FEATURED TECHNOLOGIES







#### Natural Rubber Latex

#### Latex, powder-free micro surgical gloves offer exceptional tactile sensitivity and dexterity

#### **KEY FEATURES & BENEFITS**

• Micro thin: Made 20% thinner\* for enhanced precision performance

Brown color: Reduces glare of bright surgical lights to prevent eye fatigue

• Smooth with micro-textured finish: Unique surface offers secure grip

\*When compared with Ansell's standard latex gloves

#### FEATURED TECHNOLOGIES



🛟 DERMASHIELD"

#### STANDARDS AND CERTIFICATIONS







### LATEX RANGE

### **PremierPro<sup>™</sup> Latex Ortho**

Thickness	Ortho	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Textured fi	ngers and palm
Internal Glove Surface	Polymer co DERMASHI	ated with ELD <sup>™</sup> Technology
Double Gloving Recommendation	Outer glov	e
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	6	44660
	6.5	44665
Order Number	7	44670
Order Number	7.5	44675
	8	44680
	8.5	44685
	9	44690

### **PremierPro<sup>™</sup> Latex Underglove**

Thickness	Micro	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Smooth with finish	h micro-textured
Internal Glove Surface	Coated with Technology	n HYDRASOFT™
Double Gloving Recommendation	Underglove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	5.5	44255
	6	44260
Order Number	6.5	44265
	7	44270
	7.5	44275
	8	44280
	85	44285
	0.5	44205

Latex surgical glove with an extra-thick glove formulation to withstand the rigors of surgery, especially where power tools and bone fragments are present

#### **KEY FEATURES & BENEFITS**

- **50% thicker\*:** Protection when power tools and bone fragments are present
- Heavy-textured finish: Superior grip in both dry and wet conditions
- Donning aid technology: Fast, easy donning with either dry or damp hands
- \*When compared with standard Ansell latex surgical glove alternatives

#### FEATURED TECHNOLOGIES



Latex, powder-free surgical green underglove for

Effortless donning and optimal comfort when

• Effortless dry and damp donning as well as intra-

🛟 DERMASHIELD"

SUREFIT<sup>®</sup>

**SUREFIT**<sup>®</sup>

#### STANDARDS AND CERTIFICATIONS



effortless double gloving

**KEY FEATURES & BENEFITS** 

double gloved

operative changes

FEATURED TECHNOLOGIES

STANDARDS AND CERTIFICATIONS

EN ISO 374-5

**\$** 

🔶 A.R.T

EN ISO 374-1

Type B 

+ HYDRASOFT

.

• Helps keep skin moist and supple

• Easy identification of a glove breach

### Natural Rubber Latex

Scan or <u>click here</u> for more info

Natural Rubber Latex



#### Scan or <u>click here</u> for more info



### **ENCORE®** Latex Textured

Thickness	Standard	
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology	
External Glove Surface	Textured fing	gers and palm
Internal Glove Surface	Polymer coa DERMASHIE	ted with LD™ Technology
Double Gloving Recommendation	Outer glove	
Tested for Use with Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)	
Freedom From Holes (AQL)	0.65	
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case	
	SIZE	CODE
	5.5	45755
	6	45760
Order Number	6.5	45765
	7	45770
	7.5	45775
	8	45780
	8.5	45785
	9	45790



### **ENCORE® PERRY STYLE 42**

Cuff StyleBeadedExternal Glove SurfaceSmooth with micro-textured finishInternal Glove SurfacePolymer coated	Standard	
Surface         micro-textured finish           Internal Glove         Polymer coated	Beaded	
Polymer costed		
Double Gloving Recommendation		
Tested for Use with Chemotherapy Drugs No, not tested in accordance with ASTM D6978	No, not tested in accordance with ASTM D6978	
Freedom From Holes (AQL) 0.65	0.65	
Packaging 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case		
SIZE CODE		
5.5 45355		
6 45360		
6.5 <b>45365</b>		
Order Number 7 45370		
7.5 45375		
8 45380		
8.5 45385		
9 45390		

# **KEY FEATURES & BENEFITS**

- or damp conditions
- sensitivity and protection

## 🛟 A.R.T

ERGOFORM<sup>®</sup>

2™ PF	
	Latex, powd grip for con
	KEY FEATURE
ish	Textured     or damp
	• Optimize

- sensitivity and protection

### LATEX RANGE

#### Natural Rubber Latex

#### Latex, powder-free surgical gloves ensuring secure grip for confident instrument handling

• Textured finish: Secure grip of instruments under dry

• **Optimized thickness:** Ideal balance of tactile

• Proprietary donning technology: Effortless donning and doffing with either dry or damp hands



#### FEATURED TECHNOLOGIES

DERMASHIELD<sup>®</sup>

SUREFIT





Natural Rubber Latex

#### der-free surgical gloves ensuring secure nfident instrument handling

#### RES & BENEFITS

d finish: Secure grip of instruments under dry o conditions

Optimized thickness: Ideal balance of tactile

Proprietary donning technology: Effortless donning and doffing with either dry or damp hands







# RESOURCES

PremierPro<sup>™</sup> Contract Number: PP-OR-2039/AD-OR-2039

## **GLOVE SIZE FINDER**

Getting that perfect fit is paramount, not only for comfort, but also for safety. With out glove size finder application, putting on the right size has never been easier. Our application is designed to guide users to their ideal glove size, ensuring maximum protection without compromising dexterity or ease of movement.

### Here's How It Works:






Access the application by clicking or scanning the QR Code below



Scan or <u>click here</u> for surgical gloves

Follow the simple instructions on-screen







4

### Get your glove size instantly

\*Based On The Shape Of Your Hand And The Quality Of The Scan, Your Actual Glove Size May Vary. Please Follow The Instructions Closely And We Recommend That You Try The Glove Size Suggested By The Scan.



## WHY AND HOW TO DOUBLE GLOVE?


The main causes of blood exposure by operating personnel are surgical glove punctures and tears in the operating room.<sup>1</sup> Double gloving is the recommended best practice to provide an additional layer of protection against the risk of glove microperforations and potential pathogen exposure. Be sure to follow the steps below when double gloving.

#### **Closed Self-Donning Steps**



After opening the underglove inner wrap, remove the first glove by grasping only the folded cuff edge. Lay the glove on your wrist with fingers pointing to the elbow, and with the edge of glove cuff aligned with the edge of the gown cuff. After performing surgical hand antisepsis and donning PPE as per recommended best practice, hands of the scrub team member should remain inside the gown sleeves.





AORN recommends that scrub members

Removing both pairs of gloves when a

change their gloves every 90-150 minutes.

perforation occurs is another best practice to

consider since a perforation in the outer glove

is an indication that the inner glove may be

• The edge of the gown cuff is low- towards the fingers and covering the lower palm • There is an adequate length of the glove cuff on the arm of the gown to provide grip • There is a minimal amount of folds where the glove meets the gown



#### References:

compromised as well.<sup>2</sup>

- Smoot EC, Practical precautions for avoiding sharp injuries and blood exposure. Plast Reconstr Surg., 1998; 101(2): p 528-534.
   Thomas-Copeland J. Do surgical personnel really need to double glove? AORN J. 2009;89(2): 322-328.

- Rev 2014 7:(3)



Grasp the underside of the glove cuff and gown sleeve interface together. Then, taking the front of the glove cuff, stretch over your hand, completely covering the gown cuff.



Keeping your fingers straight, slide and move fingers into the glove by pulling up the gown.



With the first pair of gloves donned and gown cuffs completely covered, attend to finer adjustments without risk of contamination, as well as inspect the gloves for perforations or tears. To prevent glove/gown separation ensure the following:



Then, open and don outer gloves by sliding fingers under the cuffs of the opposite hand gloves, and slipping over fingers and hands. Again, inspect and adjust for comfort.

3. Laine T, Aarnio P. How often does glove perforation occur in surgery? Comparison between single gloves and a double-gloving system. American journal of surgery 2001;181:564-6.

Mischke C., Verbeek J., Annika Saarto A., Lavoie MC., et al., Gloves, extra gloves or special types of gloves for preventing percutaneous exposure injuries in healthcare personnel, Cochrane Database Syst



#### Scan or <u>click here</u> for more info



Ansell North America US Tel: 800 952 9916 Email: customersolutionsus@ansell.com



Contract Number: PP-OR-2039/AD-OR-2039 TSK-30901\_Premier Surgical Brochure\_2024



PremierPro<sup>™</sup> and PremierPro<sup>™</sup> logo are trademarks of Premier Healthcare Alliance, L.P. Ansell, <sup>®</sup> and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking. © 2024 Ansell Limited. All Rights Reserved.