



## EXAM GLOVE PORTFOLIO



As a direct sourcing alternative, S2S Global benefits Premier members by removing costs from the supply chain and by providing a greater impact on product flow, pricing and quality under the PremierPro® brand.

S2S Global offers a comprehensive selection of **NITRILE**, **VINYL** and **LATEX** exam gloves to meet the many needs of healthcare facilities. Choosing the right glove is a simple process. First consider the use and then what features will be best for your staff.

**All PremierPro Exam Gloves are powder-free and ambidextrous with beaded cuffs, designed for single use only.**

Please consult the following pages for additional details specific to each glove line.



# EXAM GLOVES SPECIFICATIONS

## EXAM GLOVES COMPARISON TABLE

		Texture	Finger Thickness	Palm Thickness	ASTM 6978 <sup>†</sup>	USP 800 <sup>‡</sup>	
<b>500</b>	<b>Enhanced Chemo Tested Nitrile Exam Gloves</b>	Fingertip Textured	0.08mm	0.06mm	✓	✓	
<b>506</b>	<b>Plus Improved Thinner Technology Nitrile Exam Gloves</b>	Fingertip Textured	0.09mm	0.07mm	✓	✓	
<b>504</b>	<b>Plus Nitrile Exam Gloves</b>	Fingertip Textured	0.12mm	0.10mm	✓	✓	
<b>509</b>	<b>Extended Cuff Nitrile Exam Gloves</b>	Fingertip Textured	0.19mm	0.15mm	✓	✓	
<b>507</b> <b>508</b>	<b>Sterile Nitrile Exam Gloves</b>	Lightly Textured	0.15mm	0.13mm			
<b>403</b> <b>404</b>	<b>Vinyl Exam Gloves</b>	Smooth	0.08mm-0.10mm	0.08mm-0.10mm			
<b>406</b> <b>407</b> <b>409</b>	<b>Stretch Vinyl Exam Gloves</b>	Smooth	0.10mm-0.12mm	0.10mm			
<b>460</b>	<b>Powder-Free Textured Latex Exam Gloves</b>	Textured	0.12mm	0.12mm			

◆ Complies with standard ASTM 6978 (Assessment of Medical Exam Gloves for Resistance to Permeation by Chemotherapy Drugs).  
 ▲ Meets USP 800 Personal Protection Equipment requirement for use of Chemotherapy Gloves when handling Hazardous Drugs.

## NITRILE EXAM GLOVES

All PremierPro™ Nitrile Exam Gloves are made with Synthetic Nitrile Rubber (NBR) and feature beaded cuffs. They are powder-free, ambidextrous and designed for single use only. PremierPro™ Nitrile exam gloves are not made with natural rubber latex.




### 500 - Enhanced Chemo Tested Nitrile Exam Gloves

PremierPro™ Enhanced Chemo Tested Nitrile Exam Gloves, our newest line addition, offer a combination of tactile sensitivity, fit, strength, and value. With this glove, you have the security of knowing it has been tested for use with a total of 36 Hazardous Drugs, complies with standard ASTM 6978, meets USP 800 PPE requirements, and meets NIOSH, ONS and OSHA recommendations. This glove is ideal for wide variety of clinical tasks and procedures.

### 506 - Plus Improved Thinner Technology Nitrile Exam Gloves

Our most popular glove, the PremierPro™ Plus Improved Thinner Technology Nitrile Exam Gloves combine feel, fit, strength, comfort, tactile sensitivity and value. The PremierPro Plus line offers additional enhancements including thicker fingertips and improved donning. These gloves are tested for use with chemotherapy drugs.



### 504 - Plus Nitrile Exam Gloves

PremierPro™ Plus Nitrile Exam Gloves offer superb elasticity and ultra-comfortable fit in a glove not made with natural rubber latex. The PremierPro Plus line offers additional enhancements including thicker fingertips and improved donning. These gloves are tested for use with chemotherapy drugs.

### 509 - Extended Cuff Nitrile Exam Gloves

PremierPro™ Extended Cuff Nitrile Exam Gloves are 12 inches in length and offer superior elasticity to provide a better feel and a more comfortable fit. The extra length and extra thickness provide greater protection for hands as well as over the wrist and forearms. Designed for the health care providers' protection in mind against body fluids, and for use with Chemotherapy agents in high-risk clinical settings.



### 507/508 - Sterile Nitrile Exam Gloves

PremierPro™ Sterile Nitrile Exam Gloves offer superior elasticity to provide a better feel and a more comfortable fit. Ideal for use with sterile procedures at the patient bedside, such as endotracheal suctioning. Gloves are lightly textured which improves their ability to grip. Designed for the health care providers' protection in mind against body fluids.

VINYL EXAM GLOVES



**403/404 - Vinyl Exam Gloves**

PremierPro™ Vinyl Exam Gloves are general use gloves which are ambidextrous and not made with natural rubber latex to eliminate the risk of skin sensitivity to latex. They feature a smooth finish for greater tactile sensitivity and excellent donning characteristics. These gloves are available in **standard** and **economy** grades. Standard grade gloves are thicker than the economy grade gloves.



**406/407/409 - Stretch Vinyl Exam Gloves**

PremierPro™ Stretch Vinyl Exam Gloves are general use gloves which feature a smooth finish and a softer vinyl for more flexibility and comfort, as well as excellent donning characteristics. The smooth finish provides greater tactile sensitivity making the stretch vinyl exam glove suitable for a variety of tasks. These gloves are available in **standard**, **heavy** and **platinum** grades. Heavy stretch vinyl is a more traditional stretch vinyl glove. Platinum stretch vinyl also features an ultra soft formula for improved comfort.



LATEX EXAM GLOVES



**460 - Powder-Free Textured Latex Exam Gloves**

PremierPro™ Latex Exam Gloves are designed to provide the ideal combination of protection, tactile sensitivity and value for a variety of needs. Textured surface offers improved tactile sensitivity with enhanced grip. This glove is easy to don and is both strong and flexible.



# EXAM GLOVE PACKAGING INDEX

			X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
<b>Enhanced Chemo Tested Nitrile Exam Gloves</b>	<b>REORDER #</b>		5001	5002	5003	5004	5005
	<b>PACKAGING</b>		250/bx; 10 bxs/cs	250/bx; 10 bxs/cs	250/bx; 10 bxs/cs	250/bx; 10 bxs/cs	250/bx; 10 bxs/cs
<b>Plus Improved Thinner Technology Nitrile Exam Gloves</b>	<b>REORDER #</b>		5061	5062	5063	5064	5065
	<b>PACKAGING</b>		200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	180/bx; 10 bxs/cs
<b>Plus Nitrile Exam Gloves</b>	<b>REORDER #</b>		5041	5042	5043	5044	5045
	<b>PACKAGING</b>		200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	200/bx; 10 bxs/cs	180/bx; 10 bxs/cs
<b>Extended Cuff Nitrile Exam Gloves</b>	<b>REORDER #</b>		—	5092	5093	5094	5095
	<b>PACKAGING</b>		—	50/bx; 10 bxs/cs	50/bx; 10 bxs/cs	50/bx; 10 bxs/cs	45/bx; 10 bxs/cs
<b>Sterile Nitrile Exam Gloves</b>	<b>Single</b>	<b>REORDER #</b>	—	5072	5073	5074	5075
		<b>PACKAGING</b>	—	1 gl/pch; 100 pchs/bx; 4 bxs/cs	1 gl/pch; 100 pchs/bx; 4 bxs/cs	1 gl/pch; 100 pchs/bx; 4 bxs/cs	1 gl/pch; 100 pchs/bx; 4 bxs/cs
	<b>Pair</b>	<b>REORDER #</b>	—	5082	5083	5084	5085
		<b>PACKAGING</b>	—	1 pr/pch; 50 pchs/bx; 4 bxs/cs	1 pr/pch; 50 pchs/bx; 4 bxs/cs	1 pr/pch; 50 pchs/bx; 4 bxs/cs	1 pr/pch; 50 pchs/bx; 4 bxs/cs
<b>Vinyl Exam Gloves</b>	<b>Standard</b>	<b>REORDER #</b>	4041	4042	4043	4044	4045
		<b>PACKAGING</b>	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs
	<b>Economy</b>	<b>REORDER #</b>	—	4032	4033	4034	4035
		<b>PACKAGING</b>	—	150/bx; 10 bxs/cs	150/bx; 10 bxs/cs	150/bx; 10 bxs/cs	130/bx; 10 bxs/cs
<b>Stretch Vinyl Exam Gloves</b>	<b>Standard</b>	<b>REORDER #</b>	4071	4072	4073	4074	4075
		<b>PACKAGING</b>	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs
	<b>Heavy</b>	<b>REORDER #</b>	4061	4062	4063	4064	4065
		<b>PACKAGING</b>	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	90/bx; 10 bxs/cs
	<b>Platinum</b>	<b>REORDER #</b>	4091	4092	4093	4094	4095
		<b>PACKAGING</b>	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs
<b>Powder-Free Textured Latex Exam Gloves</b>	<b>REORDER #</b>		4601	4602	4603	4604	4605
	<b>PACKAGING</b>		100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	100/bx; 10 bxs/cs	90/bx; 10 bxs/cs

Premier Contract Numbers: PP-NS-1231, AS-NS-1231





## **DOUBLE GLOVING**

With the PremierPro™ full line of Exam Gloves and the spectrum of colors offered, you can easily double glove in high risk situations. The underglove will be clearly visible if the top layer is breached, as long as the two layers are from different lines of Exam Gloves.

Safety is a high priority at S2S Global and we are proud to provide an easy solution for challenging environments at our member hospital systems.

**Double Gloving using  
504 and 509  
PremierPro™ Exam Gloves**

**For more information about S2S Global's  
comprehensive line of Exam Gloves, please  
contact your PremierPro™ sales representative.**

**S2S GLOBAL Customer Service:**

**855.531.7699  
info@s2s-global.com**

