

HaloGUARD™
HaloGUARD™ Protective Disc


- Helps prevent contamination of the wound dressing for 7 days (168 hours)
- Chlorhexidine gluconate (CHG) Antimicrobial agent consistent with CDC and DOQI guidelines
- Broad-spectrum of antimicrobial agent
- 5 log CFU/ml reduction of 8 different target pathogens
- Pyrogen and latex free
- Sterility assurance level: 10⁻⁶
- Film adhesive: Pressure sensitive acrylic
- Softer material to help reduce skin irritation/abrasion
- Higher absorbency of foam (38% higher than market-leading product)
- Performed tests of in-vitro Cytotoxicity, irritation and skin sensitization, and chemical characterization
- Better in-vivo wound healing results

DESCRIPTION	REORDER NUMBER	PACKAGING
1 inch (2.5 cm) disc, 4.0mm center hole with radial slit	Hg01410	10 Each/Carton; 40 Cartons/Case
1 inch (2.5 cm) disc, 7.0mm center hole with radial slit	Hg01710	10 Each/Carton; 40 Cartons/Case
0.75 inch (1.9 cm) disc, 1.5mm center hole with radial slit	Hg00751510	10 Each/Carton; 40 Cartons/Case

Please refer to HaloGUARD™ Product Brochure as source of information