




PremierPro™ FLUSH SYRINGES



SAFETY

- Pre-filled with saline 0.9% NaCl for the flushing of central and peripheral vascular access devices
- Single-use, pre-filled syringes can help reduce medication errors and contamination risks associated with multi-dose vials
- Terminal sterilization for maximum sterility and safe administration
- Preservative-free & DEHP-free
- Not made with natural rubber latex



Color-coded white tip for easy visual identification

Large, accurate graduations help clinicians quickly identify the correct dose size

Smooth plunger design that gives the clinician excellent control over the syringe

CLINICIAN FRIENDLY

- Finger grip design that is comfortable and gives user more control over the syringe during flushing
- Plunger rod design to help reduce reflux
- Clear, concise, and accurate large graduations on the syringe barrel to help user deliver precise volume of solution
- Peel pouch is made of clear plastic to view the product before opening

Everything you need for catheter flushing



CONVENIENCE

- ISO- compliant syringe design to help with your facility's compliance
- Ready-to-use, pre-filled and pre-assembled; no additional components to buy
- Barcoded with QR code for bedside scanning/inventory management
- 3-year Shelf Life



Item No.	Description	Packaging
FL1003D	Pre-Filled Normal Saline Flush Syringe 3 mL Fill 0.9% Sodium Chloride Injection, USP. Sterile Fluid Pathway.	30/BX, 6 BX/CS
FL1005D	Pre-Filled Normal Saline Flush Syringe 5 mL Fill 0.9% Sodium Chloride Injection, USP. Sterile Fluid Pathway.	30/BX, 6 BX/CS
FL1010D	Pre-Filled Normal Saline Flush Syringe 10 mL Fill 0.9% Sodium Chloride Injection, USP. Sterile Fluid Pathway.	30/BX, 6 BX/CS
FL1010SF	Pre-Filled Normal Saline Flush Syringe 10 mL Fill 0.9% Sodium Chloride Injection, USP. Sterile Field Ready.	30/BX, 6 BX/CS