


## HaloGUARD™

### HaloGUARD™ Protective Disc



- Helps prevent contamination of the wound dressing for 7 days (168 hours)
- Chlorhexidine gluconate (CHG) Antimicrobial agent consistent with CDC and DOQI guidelines
- Broad-spectrum of antimicrobial agent
- 5 log CFU/ml reduction of 8 different target pathogens
- Pyrogen and latex free
- Sterility assurance level: 10-6
- Film adhesive: Pressure sensitive acrylic
- Softer material to help reduce skin irritation/abrasion
- Higher absorbency of foam (38% higher than market-leading product)
- Performed tests of in-vitro Cytotoxicity, irritation and skin sensitization, and chemical characterization
- Better in-vivo wound healing results

| DESCRIPTION   | REORDER NUMBER | PACKAGING                          |
|---|----------------|------------------------------------|
| 1 inch (2.5 cm) disc, 4.0mm center hole with radial slit    | Hg01410        | 10 Each/Carton;<br>40 Cartons/Case |
| 1 inch (2.5 cm) disc, 7.0mm center hole with radial slit    | Hg01710        | 10 Each/Carton;<br>40 Cartons/Case |
| 0.75 inch (1.9 cm) disc, 1.5mm center hole with radial slit | Hg00751510     | 10 Each/Carton;<br>40 Cartons/Case |

\*Please refer to HaloGUARD™ Product Brochure as source of information

## Pre-Filled Flush Syringes



- ISO- compliant syringe design to help with your facility's compliance
- Ready-to-use, pre-filled and pre-assembled; no additional components to buy
- Barcoded with QR code for bedside scanning/inventory management
- 3-year Shelf Life
- Pre-filled with saline 0.9% NaCl for the flushing of central and peripheral vascular access devices
- Single-use pre-filled syringes can help reduce medication errors and contamination risks associated with multi-dose vials
- Terminal sterilization for maximum sterility and safe administration
- Preservative-free & DEHP-free
- Not made with natural rubber latex

| DESCRIPTION   | REORDER NUMBER | PACKAGING     |
|---|----------------|---------------|
| Pre-Filled Normal Saline Flush Syringe<br>3 mL Fill 0.9% Sodium Chloride<br>Injection, USP. Sterile Fluid Pathway.  | FL1003D        | 30/BX, 6BX/CS |
| Pre-Filled Normal Saline Flush Syringe<br>5 mL Fill 0.9% Sodium Chloride<br>Injection, USP. Sterile Fluid Pathway.  | FL1005D        | 30/BX, 6BX/CS |
| Pre-Filled Normal Saline Flush Syringe<br>10 mL Fill 0.9% Sodium Chloride<br>Injection, USP. Sterile Fluid Pathway. | FL1010D        | 30/BX, 6BX/CS |
| Pre-Filled Normal Saline Flush Syringe<br>10 mL Fill 0.9% Sodium Chloride<br>Injection, USP. Sterile Field Ready.   | FL1010SF       | 30/BX, 6BX/CS |