S2S GLOBAL PremierPro

## **EXAM GLOVES**

## STERILE NITRILE EXAM GLOVES

PremierPro<sup>™</sup> Sterile Nitrile Exam Gloves offer superior elasticity to provide a better feel and a more comfortable fit. Ideal for use with sterile procedures at the patient bedside, such as endotracheal suctioning. Gloves are lightly textured which improves their ability to grip. Designed for the health care providers' protection in mind against body fluids. Nitrile is a synthetic material that decreases the risk of latex sensitivity in susceptible individuals.

## Features & Benefits

- Sterile
- Cuffed and walleted
- Synthetic Nitrile Rubber (NBR)
- Powder-free
- 12 Inch Beaded Cuff
- Single use only
- Ambidextrous
- Not made with natural rubber latex
- Meets USP 800 Personal Protection Equipment requirement for use of Chemotherapy Gloves when handling Hazardous Drugs



DESCRIPTION	SIZE	REORDER NUMBER	PACKAGING
Single - Sterile	Small	5072	1 gl/pch; 100 pchs/bx; 4 bxs/cs
Single - Sterile	Medium	5073	1 gl/pch; 100 pchs/bx; 4 bxs/cs
Single - Sterile	Large	5074	1 gl/pch; 100 pchs/bx; 4 bxs/cs
Single - Sterile	X-Large	5075	1 gl/pch; 100 pchs/bx; 4 bxs/cs
Pair - Sterile	Small	5082	1 pr/pch; 50 pchs/bx; 4 bxs/cs
Pair - Sterile	Medium	5083	1 pr/pch; 50 pchs/bx; 4 bxs/cs
Pair - Sterile	Large	5084	1 pr/pch; 50 pchs/bx; 4 bxs/cs
Pair - Sterile	X-Large	5085	1 pr/pch; 50 pchs/bx; 4 bxs/cs
Premier Contract Numbers: PD-NS-1231 AD-NS-1231 SD-NS-1869			

Premier Contract Numbers: PP-NS-1231, AD-NS-1231, SP-NS-1869

AVERAGE NORMALIZED BREAKTHROUGH TIME IN MINUTES
>240
>240
>240
>240
>240
>240
>240
>240
>240
>240

Warning: Do not use with Carmustine or Thiotepa

Physical Properties (Before and After Aging) Tensile Strength: 14 Mpa minimum Elongation: 500% minimum Length: 290mm Finger Thickness: 0.14mm Palm Thickness: 0.09mm